

TRIP6 Antibody



PACO12929

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

TRIP6

Uniprot

Q15654

Synonyms:

thyroid hormone receptor interactor 6; TRIP6; MGC10556; MGC10558; MGC29959; MGC3837; MGC4423; OIP1; ZRP-1 ;

Immunogen:

Human TRIP6.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A