

# TROAP Antibody



PACO12934

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TROAP

**Uniprot**

Q12815

**Synonyms:**

trophinin associated protein (tastin); TROAP; TASTIN ;

**Immunogen:**

Human TROAP.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A