

TRPA1 Antibody



PACO12936

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TRPA1

Uniprot

O75762

Synonyms:

transient receptor potential cation channel, subfamily A, member 1; TRPA1; ANKTM1 ;

Immunogen:

Human TRPA1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A