

TRPC1 Antibody



PACO12937

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TRPC1

Uniprot

P48995

Synonyms:

transient receptor potential cation channel, subfamily C, member 1; TRPC1; HTRP-1; MGC133334; MGC133335; TRP1 ;

Immunogen:

Human TRPC1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A