

# TRPC4AP Antibody



PACO12939

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TRPC4AP

**Uniprot**

Q8TEL6

**Synonyms:**

transient receptor potential cation channel, subfamily C, member 4 associated protein; TRPC4AP; C20orf188; TRRP4AP; TRUSS ;

**Immunogen:**

Human TRPC4AP.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A