

TSC22D1 Antibody



PACO12947

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

TSC22D1

Uniprot

Q15714

Synonyms:

16kd;

Immunogen:

Human TSC22D1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.320° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A