

TSEN2 Antibody



PACO12950

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TSEN2

Uniprot

Q8NCE0

Synonyms:

tRNA splicing endonuclease 2 homolog (S. cerevisiae); TSEN2; MGC2776; MGC4440; PCH2B; SEN2; SEN2L ;

Immunogen:

Human TSEN2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A