

# TSHZ1 Antibody



PACO12955

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TSHZ1

**Uniprot**

Q6ZSZ6

**Synonyms:**

teashirt zinc finger homeobox 1; TSHZ1; NY-CO-33; SDCCAG33; TSH1 ;

**Immunogen:**

Human TSHZ1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A