

# TSNAX Antibody



PACO12957

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:****Protein Background:****Gene ID:**

TSNAX

**Uniprot**

Q99598

**Synonyms:**

translin-associated factor X; TSNAX; TRAX ;

**Immunogen:**

Human TSNAX.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A