

TSPAN12 Antibody



PACO12958

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TSPAN12

Uniprot

O95859

Synonyms:

tetraspanin 12; TSPAN12; NET-2; TM4SF12 ;

Immunogen:

Human TSPAN12.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A