

TSPAN13 Antibody



PACO12959

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TSPAN13

Uniprot

O95857

Synonyms:

tetraspanin 13; TSPAN13; FLJ22934; NET-6; TM4SF13 ; Tspan-13

Immunogen:

Human TSPAN13.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A