

# TSPAN5 Antibody



PACO12962

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TSPAN5

**Uniprot**

P62079

**Synonyms:**

tetraspanin 5; TSPAN5; NET-4; TM4SF9; TSPAN-5 ;

**Immunogen:**

Human TSPAN5.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A