

# TSPYL4 Antibody



PACO12965

---

## Product Information

**Size:**

50ul

**Reactivity:**

Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TSPYL4

**Uniprot**

Q9UJ04

**Synonyms:**

TSPY-like 4; TSPYL4; KIAA0721; dJ486I3.2 ;

**Immunogen:**

Human TSPYL4.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A