

# TSR2 Antibody



PACO12967

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TSR2

**Uniprot**

Q969E8

**Synonyms:**

TSR2, 20S rRNA accumulation, homolog (S. cerevisiae); TSR2; DT1P1A10; MGC20451; WGG1 ;

**Immunogen:**

Human TSR2.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A