

TSTA3 Antibody



PACO12969

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TSTA3

Uniprot

Q13630

Synonyms:

tissue specific transplantation antigen P35B; TSTA3; FX; P35B; SDR4E1 ;

Immunogen:

Human TSTA3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A