

TTC1 Antibody



PACO12971

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

TTC1

Uniprot

Q99614

Synonyms:

tetratricopeptide repeat domain 1; TTC1; FLJ46404; TPR1 ;

Immunogen:

Human TTC1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A