

TUBA1B Antibody



PACO12982

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

TUBA1B

Uniprot

P68363

Synonyms:

tubulin, alpha 1b; TUBA1B; K-ALPHA-1 ; alpha-tubulin, a-tubulin

Immunogen:

Human TUBA1B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; Avoid freeze / thaw cycles.

Product Images

N/A

N/A