

TUBGCP4 Antibody



PACO12989

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TUBGCP4

Uniprot

Q9UGJ1

Synonyms:

tubulin, gamma complex associated protein 4; TUBGCP4;76P; FLJ14797; GCP4 ;

Immunogen:

Human TUBGCP4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A