

TULP4 Antibody



PACO12993

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

TULP4

Uniprot

Q9NRJ4

Synonyms:

tubby like protein 4; TULP4; KIAA1397; TUSP ;

Immunogen:

Human TULP4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A