

PACO12996

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

TXK

**Uniprot**

P42681

**Synonyms:**

TXK tyrosine kinase; TXK; BTKL; MGC22473; PSCTK5; PTK4; RLK; TKL ;

**Immunogen:**

Human TXK.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A