

PACO12997

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

TXN

Uniprot

P10599

Synonyms:

thioredoxin; TXN; DKFZp686B1993; MGC61975; TRX; TRX1 ;

Immunogen:

Human TXN.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A