

TXNDC17 Antibody



PACO13000

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

TXNDC17

Uniprot

Q9BRA2

Synonyms:

thioredoxin domain containing 17; TXNDC17; MGC14353; TRP14; TXNL5 ;

Immunogen:

Human TXNDC17.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A