

NME8 Antibody



PACO13001

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

NME8

Uniprot

Q8N427

Synonyms:

thioredoxin domain containing 3 (spermatozoa); TXNDC3; CILD6; NME8; SPTRX2 ;

Immunogen:

Human TXNDC3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A