

PACO13007

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

TYMS

Uniprot

P04818

Synonyms:

thymidylate synthetase; TYMS; HsT422; MGC88736; TMS; TS; TSase ;

Immunogen:

Human TYMS.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A