

UBA2 Antibody



PACO13010

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

UBA2

Uniprot

Q9UBT2

Synonyms:

ubiquitin-like modifier activating enzyme 2; UBA2; ARX; FLJ13058; HRIHFB2115; SAE2 ;

Immunogen:

Human UBA2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A