

# UBA52 Antibody



PACO13012

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

UBA52

**Uniprot**

P62987

**Synonyms:**

ubiquitin A-52 residue ribosomal protein fusion product 1; UBA52; CEP52; HUBCEP52; MGC126879; MGC126881; MGC57125; RPL40 ;

**Immunogen:**

Human UBA52.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A