

UBA7 Antibody



PACO13014

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

UBA7

Uniprot

P41226

Synonyms:

ubiquitin-like modifier activating enzyme 7; UBA7; D8; MGC12713; UBA1B; UBE1L; UBE2
;

Immunogen:

Human UBA7.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A