

UBE2E2 Antibody



PACO13026

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

UBE2E2

Uniprot

Q96LR5

Synonyms:

ubiquitin-conjugating enzyme E2E 2 (UBC4/5 homolog, yeast); UBE2E2; FLJ25157; UBCH8 ;

Immunogen:

Human UBE2E2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A