

UBE2G2 Antibody



PACO13030

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

UBE2G2

Uniprot

P60604

Synonyms:

ubiquitin-conjugating enzyme E2G 2 (UBC7 homolog, yeast); UBE2G2; UBC7 ;

Immunogen:

Human UBE2G2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A