

UBE2L3 Antibody



PACO13033

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

UBE2L3

Uniprot

P68036

Synonyms:

ubiquitin-conjugating enzyme E2L 3; UBE2L3; E2-F1; L-UBC; UBCH7; UbcM4 ;

Immunogen:

Human UBE2L3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A