

UBE2N Antibody



PACO13035

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

UBE2N

Uniprot

P61088

Synonyms:

ubiquitin-conjugating enzyme E2N (UBC13 homolog, yeast); UBE2N; MGC131857; MGC8489; UBC13; UbcH-ben ;

Immunogen:

Human UBE2N.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A