

UBE2O Antibody



PACO13036

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

UBE2O

Uniprot

Q9C0C9

Synonyms:

ubiquitin-conjugating enzyme E2O; UBE2O; E2-230K; FLJ12878; KIAA1734 ;

Immunogen:

Human UBE2O.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A