

UBE2U Antibody



PACO13040

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

UBE2U

Uniprot

Q5VWX9

Synonyms:

ubiquitin-conjugating enzyme E2U (putative); UBE2U; MGC35130; RP4-636O23.1 ;

Immunogen:

Human UBE2U.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A