

UBE2V2 Antibody



PACO13042

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

UBE2V2

Uniprot

Q15819

Synonyms:

ubiquitin-conjugating enzyme E2 variant 2; UBE2V2; DDVIT1; DDV'it-1; EDAF-1; EDPF-1; EDPF1; MMS2; UEV-2; UEV2 ;

Immunogen:

Human UBE2V2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A