

UBE2W Antibody



PACO13043

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

UBE2W

Uniprot

Q96B02

Synonyms:

ubiquitin-conjugating enzyme E2W (putative); UBE2W; FLJ11011; hUBC-16 ;

Immunogen:

Human UBE2W.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A