

UBL4B Antibody



PACO13050

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

UBL4B

Uniprot

Q8N7F7

Synonyms:

ubiquitin-like 4B; UBL4B; FLJ25690 ;

Immunogen:

Human UBL4B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A