

# UBL7 Antibody



PACO13051

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

UBL7

**Uniprot**

Q96S82

**Synonyms:**

ubiquitin-like 7 (bone marrow stromal cell-derived); UBL7; BMSC-UbP; MGC14421; TCBA1 ;

**Immunogen:**

Human UBL7.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A