

UBLCP1 Antibody



PACO13052

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

UBLCP1

Uniprot

Q8WVY7

Synonyms:

ubiquitin-like domain containing CTD phosphatase 1; UBLCP1; CPUB1; FLJ25267; MGC10067 ;

Immunogen:

Human UBLCP1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A