

UBQLN3 Antibody



PACO13055

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

UBQLN3

Uniprot

Q9H347

Synonyms:

ubiquilin 3; UBQLN3; TUP-1 ;

Immunogen:

Human UBQLN3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A