

PACO13056

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

UBR2

Uniprot

Q8IWW8

Synonyms:

ubiquitin protein ligase E3 component n-recognin 2; UBR2; C6orf133; DKFZp686C08114; KIAA0349; MGC71112; RP3-392M17.3; bA49A4.1; dJ242G1.1; dJ392M17.3 ;

Immunogen:

Human UBR2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A