

UBXN1 Antibody



PACO13059

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

UBXN1

Uniprot

Q04323

Synonyms:

UBX domain protein 1; UBXN1;2B28; UBXD10 ;

Immunogen:

Human UBXN1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A