

UBXN4 Antibody



PACO13061

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IF

Recommended dilutions:

Protein Background:

Gene ID:

UBXN4

Uniprot

Q92575

Synonyms:

UBX domain protein 4; UBXN4; FLJ23318; KIAA0242; KIAA2042; UBXD2; UBXDC1; erasin
;

Immunogen:

Human UBXN4.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A