

UCHL3 Antibody



PACO13064

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

UCHL3

Uniprot

P15374

Synonyms:

ubiquitin carboxyl-terminal esterase L3 (ubiquitin thiolesterase); UCHL3; UCH-L3 ;

Immunogen:

Human UCHL3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A