

UCHL5 Antibody



PACO13065

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

UCHL5

Uniprot

Q9Y5K5

Synonyms:

ubiquitin carboxyl-terminal hydrolase L5; UCHL5; CGI-70; INO80R; UCH37 ;

Immunogen:

Human UCHL5.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A