

UGGT2 Antibody



PACO13077

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

UGGT2

Uniprot

Q9NYU1

Synonyms:

UDP-glucose ceramide glucosyltransferase-like 2; UGCGL2; FLJ10873; FLJ11485; HUGT2; MGC117360; MGC150689; MGC87276 ;

Immunogen:

Human UGCGL2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A