

UGT2B7 Antibody



PACO13082

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

UGT2B7

Uniprot

P16662

Synonyms:

UDP glucuronosyltransferase 2 family, polypeptide B7; UGT2B7; UGT2B9 ;

Immunogen:

Human UGT2B7.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A