

# UHRF1 Antibody



PACO13087

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

UHRF1

**Uniprot**

Q96T88

**Synonyms:**

ubiquitin-like with PHD and ring finger domains 1; UHRF1; FLJ21925; ICBP90; MGC138707; Np95; RNF106; hNP95 ;

**Immunogen:**

Human UHRF1.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A