

UNC13D Antibody



PACO13093

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

UNC13D

Uniprot

Q70J99

Synonyms:

unc-13 homolog D (C. elegans); UNC13D; FHL3; HLH3; HPLH3; Munc13-4 ;

Immunogen:

Human UNC13D.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A