

# USP13 Antibody



PACO13118

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

USP13

**Uniprot**

Q92995

**Synonyms:**

ubiquitin specific peptidase 13 (isopeptidase T-3); USP13; ISOT3; IsoT-3 ;

**Immunogen:**

Human USP13.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A