

# USP14 Antibody



PACO13119

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

USP14

**Uniprot**

P54578

**Synonyms:**

ubiquitin specific peptidase 14 (tRNA-guanine transglycosylase); USP14; TGT ;

**Immunogen:**

Human USP14.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A